



© AUTOMATIC REFRACTOR © AUTOMATIC KERATOMETER © NON-CONTACT TONOMETER



 **MARCO**

THE LEADER IN VISION DIAGNOSTICS®

**MARCO <sup>3</sup> WORLD'S FIRST SYSTEM COMBINING TONOMETRY WITH AUTOREFRACTION**

**New Non-Contact Tonometer, Automatic Refractor, and Automatic Keratometer <sup>3</sup>in1**

Combining the measurement of refractive power, corneal curvature and intraocular pressure (IOP), the device provides fast and accurate ophthalmic patient care.

**SPACE SAVING**

**Compact & Efficient**

The M3 provides auto refraction, keratometry, and tonometry, while saving space and offering greater efficiency, as well as patient convenience. It eliminates the need for multiple instruments. Also, the M3 interfaces with EMR.

**FAST AND ACCURATE**

**High-Speed & Reliable Measurement**


With the latest advanced technologies, the M3 provides highly accurate, reliable measurements with significantly reduced measurement time.

**FULLY AUTOMATIC**

**EyeTracking System**

The M3 is completely automatic, using three axes: alignment, tracking, and focusing.

When the alignment is in the proper position, measurement begins automatically.



**USER FRIENDLY**

**Adjustable Monitor**

The clear 5-inch color LCD monitor with tilting function offers easy operation. If the operator needs to stand to lift the patient's eyelid, the monitor will still be able to be viewed.



TILTABLE LCD

**New Easy-Load Paper**

Printer provides fast and automatic paper loading, cutting and detachment.



EASY-LOAD PAPER

Simply open door, drop in paper, and go. No spool or feed to put paper through.

**Motorized Chinrest**

Up/Down buttons are used to adjust the motorized chinrest to the correct height for patient measurement.



ADJUSTABLE MOTORIZED CHINREST

**BALLOON TARGET**

The balloon picture is an infinity target. Automatic fogging minimizes accommodation and maintains fog through all measurement readings.



BALLOON TARGET



TARGET WITH FOG



COLOR ALIGNMENT INDICATORS



AUTO MEASUREMENT



RING DISPLAY



TONOMETRY READINGS

**PATIENT FRIENDLY**

**Gentle Air**

The new M3 has an advanced APC (Automatic Pressure Control) function which provides a softer puff of air for IOP measurement and increased patient comfort.

**Patient Detection Sensor**

When a patient is in front of the unit, the chinrest will not move during calibration and initialization.

**ADDITIONAL MEASUREMENTS**

**SLD vs. LED Technology.**

Competitors' utilization of LED technology to measure the refractive power of the optical system presents disadvantages including poor penetration of the LED through media opacities, limited data and accuracy of patients measured.

**What's the difference?**

**SLD vs. LED Technology**  
Comparison of Images on CCD Through Cataract SLD vs. Cataract LED

SLD      LED

**IN-HOUSE TRIAL DATA (MODEL EYE)**

**Cataract Images From the M3 of the Human Eye**

The Marco M3 adopts the advanced Pupil Zone Imaging Method for refraction measurement, which analyzes a wider area (max.  $\phi 4\text{mm}$ ) to obtain more reliable and realistic data that is closer to manifest refraction.

**Examples of Myopia, Hyperopia, and Astigmatism**

-6D      +6D      6D CYL

**NON-CONTACT TONOMETRY**

The M3 non-contact tonometry requires no sterilization. The automatic NCT measurements are operator independent, repeatable, and reproducible. Training is quick and easy.



**PRINTOUT**

The Automatic Refractor, Automatic Keratometer, and Non-Contact Tonometer data are automatically separated when printing. The printer comes equipped with an automatic paper cutter, eliminating torn printouts.

**SAMPLE PRINTOUT**

```

-----0002-----
ID  12345678901234567890
NAME                               M/F
JULY/28/2008                       4:10 PM
VD=12.00mm
WD=35cm
<R>  S      C      A
    - 1.75 -0.50 173 8
    - 1.25 -1.00 177 8
    - 1.25 -1.00   5 8
<-  1.25 -1.00 177>
<-  2.00 SE                >
    
```

**ENHANCED COMPREHENSIVE DIAGNOSTICS**

- Vertex Distance
- Near Working Distance
- Automatic Refraction
- Confidence Index
- SE Value

**Eyeprint**

```

TL  - 1.25 -1.00 177
CL  - 1.25 -1.00 177
    - 1.75 SE
    
```

- Thin Lens Data
- CL Conversion Data
- Spherical Equivalent

```

<R1  7.98 42.25 174>
<R2  7.65 44.00  84>
<AVE  7.82 43.25   >
<CYL  -1.75 174>
CS  12.5      PS  5.5
      (LAMP=ON)
    
```

- Automatic Keratometry
  - Sagittal K readings with eccentricity values.
- Pupil Size
- Corneal Size
  - Measurements are taken from limbus to limbus and recorded on a printout.

PD 63      N 59

IOP (mmHg)	
[R]	[L]
13	13
13	13
13	13
-----	
Avg. 13.0	13.0

- Near PD
- Pupillary Distance
- Non-Contact Tonometry

Actual printout includes data for both right and left eyes.





## SPECIFICATIONS

### ◎ AUTO REFRACTOR/KERATOMETER

Measurable Range	
Sphere	-30.00D to +25.00D (V.D. =12mm),(0.01/0.12/0.25D increments)
Cylinder	0D to ±12D (0.01/0.12/0.25D increments)
Axis	0° to 180° (1°/5° increments)
Measurable Minimum	
Pupil Diameter	2mm
Chart	Scenery chart (balloon target)
Radius Curvature	5.00 to 13.00mm (0.01mm increments)
Refractive power	25.96D to 67.50D (n=1.3375),(0.01/0.12/0.25D increments)
Astigmatism	0D to ±12.00D (0.01/0.12/0.25D increments)
Axis	0° to 180° (1°/5° increments)
Auto Tracking & Auto Shooting	X-Y-Z direction, Auto Shooting
PD Measurable Range	30.00mm to 85.00mm (indication increments: 1mm)
Corneal Size (CS)	8.1mm-14.6mm (indication increments: 0.1mm)
Pupil Measurable Range	0.8mm-12.1mm (indication increments: 0.1mm)

### ◎ NON-CONTACT TONOMETER

Measurement Range	1mmHg to 60mmHg
Measurement Range	
Setting	APC40, APC60 (APC=Automatic Puff Control)
Working Distance	11.0mm
Eye Fixation	Inner Fixation Light
Auto Tracking & Auto Shooting	X-Y-Z direction, Auto Shooting

### ◎ GENERAL INFORMATION

Monitor	Tiltable 5.7 inch color LCD
Printer	Thermal line printer with automatic paper cutter
Power Supply	AC100-240 V±10%
Power Consumption	100VA
Dimensions & Weight	10.23"(W) x 18.93"(D) x 19.9"(H) / 50.7 lbs.
Standard Accessories	Spare printer paper, Chinrest paper, Fixing pins, Power cord, Dust cover, Model eye.

The TONOREFII (M3) is manufactured by Nidek for Marco.



**EMR**

[www.mercoframes.com](http://www.mercoframes.com)



MERCOFRAMES