

# ARK-SERIES

AUTOMATIC REFRACTOR © AUTOMATIC KERATOMETER



 **MARCO**


THE LEADER IN VISION DIAGNOSTICS®

## ARK-560A UNIQUE FEATURES

### ◎ ARK-560A VIRTUAL VISION COMPARISON

The recall function allows the operator to demonstrate a comparison between the patients' unaided or glasses\* visual acuity with a subjectively corrected visual acuity.

\*Capable of data transfer with a Marco Auto Lensmeter.  
Spherical lenses only for the ARK-560A.



The VA Measurement provides the doctor with a quick assessment of how the AR measurement will affect a patient's quality of vision.

### ◎ STANDARD ARK FEATURES

- ◎ Pupil Zone Imaging Method
- ◎ SLD (Super Luminescent Diode)
- ◎ Wide Measurement Range: -30D to +25D
- ◎ Small Pupil Size Measurement
- ◎ 3D Auto Tracking and Capture
- ◎ Corneal Size Measurements
- ◎ 5.7" Tiltable Color LCD Monitor
- ◎ Easy Loading and Auto Paper Cutter


### ◎ APPLICATIONS:

- ◎ Post-op surgical follow up for cataract and LASIK patients
- ◎ Over refraction for contact lens patients
- ◎ Comparison eye screening when integrated with a Marco lensmeter

### ◎ SUBJECTIVE REFINEMENT VALUE


This feature is best used when a traditional refraction may not be able to be performed. Examples:

- ◎ Nursing Home Evaluations
- ◎ Remote Screening locations
- ◎ Vision Fairs
- ◎ School Screenings




AVAILABLE ONLY ON ARK-560A

BCVA




AVAILABLE ON ARK-530A/560A

CORNEAL SIZE



AVAILABLE ON ARK-530A/560A

PUPIL SIZE



AVAILABLE ON ARK-530A/560A

SAGITTAL K DISPLAYED

## MARCO ARK-530A/560A

**The Accuracy You Expect From Marco, With Speed, Efficiency, and Innovation.**

### © ACCURATE


**Incorporating SLD (Super Luminescent Diode).** The ARK-SERIES offers superior Zonal Ring-Image Method Technology using SLD. The first instrument of its kind to accurately and quickly measure patients with cataracts, corneal opacities, IOLs, and post LASIK.

**Sagittal K readings with eccentricity values.** Eccentricity measurements are becoming more widely utilized in CRT, pre- and post-LASIK evaluations, and contact lens fittings.

**Confidence index for each measurement.** All measurements are rated with a confidence index. Even when in the IOL mode, the measurements have a confidence index which helps to further clarify the reliability of the data obtained.

### © FAST

**EyeTracking System.** This totally automatic system provides 3 axes: alignment, tracking and focusing.



### Automatic fogging and high-speed measurements.

The ARK-SERIES's automatic fogging minimizes accommodation and maintains fog throughout all measurement readings. This saves valuable time and is ideal for children and patients who find it difficult to fixate.

Start	Fogging	Measurement	Start	Fogging	Measurement	Start	Fogging	Measurement
-------	---------	-------------	-------	---------	-------------	-------	---------	-------------


Conventional auto refractors employ a repeating sequence of fogging and measurement.

Start	Fogging	Measurement	Measurement	Measurement
-------	---------	-------------	-------------	-------------

The ARK-SERIES maintains fogging throughout multiple measurements.



COLOR ALIGNMENT INDICATORS



AUTO MEASUREMENT



RING DISPLAY





**Compare auto refraction with or without ADD power.** The ARK-SERIES includes this valuable feature to help determine whether the patient needs an exam for near vision. An easy parameter setting allows the working distance to be set between 35 and 70cm (5cm increments) or 14 and 28 inches (2 inch increments).

© **ADDITIONAL MEASUREMENTS**

**SLD vs. LED Technology.**




**SLD vs. LED Technology**  
Comparison of Images on CCD Through Cataract SLD vs. Cataract LED






SLD
LED




**IN-HOUSE TRIAL DATA (MODEL EYE)**

**Cataract Images From the ARK-SERIES of the Human Eye**

**Examples of Myopia, Hyperopia, and Astigmatism**

-6D
+6D
6D CYL

SLD technology measures patients with small pupils and media opacities. Competitors' utilization of LED technology to measure the refractive power of the optical system presents disadvantages including poor penetration of the LED through media opacities, limited data and diminished accuracy of patients measured.

© **PATIENT FRIENDLY**

The Marco ARK-SERIES adopts the advanced Pupil Zone Imaging Method for refraction measurement, which analyzes a wider area (max. 4mm) to obtain more reliable and realistic data that is closer to manifest refraction.

© **USER FRIENDLY**

**Adjustable Monitor**

The clear 5-inch color LCD monitor with tilting function offers easy operation. If the operator needs to stand to lift the patient's eyelid, the monitor can still be viewed.



**TILTABLE LCD**

**New Easy-Load Paper**

Printer provides fast and automatic paper loading, cutting and detachment. Simply open door, drop in paper, and go. No spool or paper feed needed.



**EASY-LOAD PAPER**

**Motorized Chinrest**

Up/Down buttons are used to adjust the motorized chinrest to the correct height for patient measurement.



**ADJUSTABLE MOTORIZED CHINREST**

© **INTELLIGENT**

**Virtual vision comparison: near and distance vision.**

Now, with the vision comparison button, the Marco LM-SERIES Automatic Lensmeter can be attached to the ARK-SERIES allowing you to compare the data from the lensmeter to the auto-refraction data. This allows the operator to ask the patient if the balloon is clear or not, providing important subjective information.



**BALLOON TARGET    FOGGED TARGET**

© MARCO CONNECT




-----0002-----

NAME M/F  
FEB/15/2009 4:10 PM

VD=13.75mm ①

WD=5cm ②

<R>	S	C	A	
	- 1.75	-0.50	173	8
	- 1.25	-1.00	177	8
	- 1.25	-1.00	5	8
<-	1.25	-1.00	177>	③
<-	2.00	SE		> ④



CL - 1.25 -1.00 177 ⑤

<R1 7.98 42.25 174>

<R2 7.65 44.00 84>

<AVE 7.82 43.25 >

<CYL -1.75 174>

<SAGITTAL> FIX ANGLE=25

SUP.	INF.	TEM.	MAS.	⑫
A 7.82	7.88	7.69	7.57	

<TOPOMETRY OF CORNEA>

SUP.	INF.	TEM.	MAS.	
S 7.56	7.64	7.89	7.81	

e	+0.21	+0.36	+0.44	+0.17	
---	-------	-------	-------	-------	--

eh	=	+0.34	Rh	=	7.84	⑬
----	---	-------	----	---	------	---

ew	=	+0.30	Rw	=	7.59
----	---	-------	----	---	------


E	=	+0.32	Ro	=	7.67
---	---	-------	----	---	------

ASTc	=	-1.38	ARa	=	+0.24
------	---	-------	-----	---	-------

ASTp	=	-1.47	dAST	=	+0.09
------	---	-------	------	---	-------

⑦ CS 12.5 PS 5.5 ⑧  
(LAMP=ON)

<L>	S	C	A	
	- 0.50	-0.75	2	9
	- 0.50	-0.75	6	9
	- 0.50	-0.75	2	8
<-	0.50	-0.75	2>	③
<-	1.00	SE		> ④




CL - 0.50 -0.75 2 ⑤

⑦ CS 12.0 PS 6.0 ⑧

⑨ PD 68 ⑩ N 64

NIDEK ARK-530A

- 1 Vertex Distance
- 2 Near Working Distance
- 3 AR Median Value
- 4 SE Value
- 5 Contact Lens Data
- 6 KM Median Value
- 7 CS Measurement
- 8 PS Measurement
- 9 Distance PD
- 10 Near PD
- 11 Eyeprint
- 12 Sagittal Radius
- 13 Eccentricity Values

 = A sample printout when the sagittal radius parameter is used.

Measurement values are printed out using a high-speed printer. Trial lens and contact lens data can be included, as can eyeprint illustrations to aid in explaining myopia, hyperopia or astigmatism to the patient. Eyeprint-only printout is available.





## ARK-SERIES SPECIFICATIONS

### ◎ AUTO REFRACTOR

Measurable Range	
<b>Sphere</b>	-30.00D to +25.00D (V.D. =12.00mm)
<b>Cylinder</b>	0D to ±12D (0.01/0.12/0.25D increments)
<b>Axis</b>	0° to 180° (1°/5° increments)
Measurable Minimum Pupil Diameter	2.00mm
Chart	Scenery chart (balloon target)

### ◎ AUTO KERATOMETER

Measurable Range	
<b>Radius Curvature</b>	5.00mm to 13.00mm (0.01mm increments)
<b>Refractive power</b>	25.96D to 67.50D (n=1.3375), (0.01/0.12/0.25D increments)
<b>Astigmatism</b>	0D to ±12.00D (0.01/0.12/0.25D increments)
<b>Axis</b>	0° to 180° (1°/5° increments)
Ordinary Measurement Area	3.30mm (R=7.70mm)
Peripheral Measurement Area	6.00mm (R=7.7mm)
Sagittal Radius Measurement	25° each from the center (Superior side, Inferior side, Temporal side, Nasal side)
VA Measurement	
Measurable Range	ARK-560A: Less than 20/200, 20/200, 20/80, 20/60, 20/50, 20/40, 20/30, 20/25, 20/20 (spherical refinement)
Correction Range	
<b>Sphere</b>	-20.00D to +20.00D (VD=12.00mm) (0.25D increments)
<b>Cylinder</b>	0.00D to +8.00D (0.25D increments)
<b>Axis</b>	0° to 180° (1°/5° increments)
PD Measurable Range	30.00mm to 85.00mm (0.10mm increments) (Near point PD: 28.00 to 80.00mm at WD=40cm)
Corneal Size Measurable Range	10.00mm to 14.00mm (0.10mm increments)
Pupil Size Measurable Range	1.00mm to 10.00mm (0.10mm increments)
Auto Tracking & Auto Shooting	X-Y-Z direction, Auto Shooting
Vision Comparison	Corrected vision with spherical and cylindrical lenses

### ◎ GENERAL INFORMATION

Monitor	Tiltable 5.7 inch color LCD
Printer	Built-in thermal type line printer (Easy loading and auto paper cutter)
Interface	RS-232C (IN/OUT), USB, IC Card Reader system* *Card is optional
Power Supply	AC100-240 V±10%, 50/60 Hz
Power Consumption	100VA
Dimensions & Weight	10.23"(W) x 18.93"(D) x 17.9"(H) / 44.10 lbs.
Standard Accessories	Spare printer paper, Chinrest paper, Fixing pins, Power cord, Dust cover, Model eye.



The ARK-530A/560A are manufactured by Nidek.

