Health o meter[®]
P r o f e s s i o n a l

PROPLUS[®] Display Module

**USER INSTRUCTIONS
INSTRUCCIONES DE USO
MODE D'EMPLOI**

**ISO13485:2003
CERTIFIED**

User Manual
PROPLUS® Display Module

TABLE OF CONTENTS

ENGLISH	3
ESPAÑOL	28
FRANÇAIS	53


User Manual
PROPLUS[®] Display Module

Thank you for your purchase of this Health o meter[®] Professional product. Please read this manual carefully and keep it handy for ready reference or training support.

TABLE OF CONTENTS

KEYPAD DESCRIPTION	4
CAUTION AND WARNING	5
SPECIFICATIONS	5
CERTIFICATIONS / CONNECTIVITY	6
BATTERY INSTALLATION	7
AC ADAPTER USE	7
QUICK START INSTRUCTIONS	7
OPERATION INSTRUCTIONS	9
MAINTENANCE	20
EXPLODED DIAGRAM	20
DISPLAY MODULE PARTS LIST	21
CALIBRATION PROCEDURE	22
PC COMMUNICATION PROTOCOL	22
WARRANTY AND SERVICE SUPPORT	28

KEYPAD DESCRIPTION

4021, 4025, 4026, 4650

CAUTION AND WARNING

To prevent injury and damage to your display module, please follow these instructions very carefully.

- Do not use in the presence of flammable materials.
- Operating at other voltages and frequencies than specified could damage the equipment.
- If the —LOW BAT” indicator activates, for accurate weighing replace the batteries or connect the scale to an AC power source as soon as possible.
- It is intended that this equipment be used with assistance of a healthcare worker.

SPECIFICATIONS

Health o meter® Professional ProPlus® Electronic Scales are highly sophisticated microprocessor technology. Each precision instrument is designed to provide accurate, reliable and repeatable weight measurements and features that make the weighing process simple, fast and convenient.

The scales are set up to use motion sensing technology, to determine the actual weight of a moving patient. The display settings may be changed to measure in live weight. See page 16 for instructions on changing the display settings.

The weight can be displayed in pounds (decimal; fractions of a pound; lb/oz) or in kilograms.

The display unit can be operated using its AC adapter (included) or by 6-D cell batteries (not included).

DISPLAY SPECIFICATIONS

<u>Capacity and Resolution</u>	
4021, 4025, 4026, 4261, 4650, and 4700	1,000 lb x 0.2 lb / 1/4 lb / 4 oz (454 kg x 0.1 kg)
4024	800 lb x 0.2 lb / 1/4 lb / 4 oz (360 kg x 0.1 kg)
4011	700 lb x 0.2 lb / 1/4 lb / 4oz (310 kg x 0.1 kg)
4020	400 lb x 0.2 lb / 1/4 lb / 4 oz (181 kb x 0.1 kg)
<u>Power Requirements</u>	Adapter model ADPT 31* (Included) 120-240VAC, 50-60 HZ — 9VDC (6) D Cell Batteries (not Included)
<u>Environmental</u>	
Operating Temperatures:	50°F to 95°F
Storage Temperatures:	32°F to 122°F
Humidity (Maximum):	85% RH

* Use only a Health o meter® Professional provided power supply.

CERTIFICATIONS / CONNECTIVITY

CERTIFICATION DESCRIPTIONS

**ISO13485:2003
CERTIFIED**

Pelstar, LLC has been officially certified as a manufacturer of Health o meter® Professional medical devices to be in compliance with stringent legal requirements. The Pelstar, LLC quality assurance system covers the development, production, sales and service of medical scales and measuring systems.

PROPLUS® SERIES SCALE CONNECTIVITY INFORMATION

This scale is designed to connect to a computer, printer, monitor or other devices via USB or serial port, to allow reliable transmission of weight, height and BMI from the scale. This connectivity helps close a source of potential medical errors caused by incorrectly copying and then recording weight information taken. This scale is EMR-EHR compatible. For more information regarding how to connect this scale to devices please refer to manual or contact Customer service through one of the numbers listed below.


CUSTOMER SERVICE INFORMATION

Our fully staffed, in-house customer service team is equipped to answer any questions or concerns that may arise.

BATTERY INSTALLATION / REPLACEMENT


1. Unplug the ProPlus® display from all power sources.
2. Using a Phillips screwdriver remove the battery cover from the display (A).
3. Disconnect battery holder cable connector from the scale-battery connector (B).
4. Carefully remove the battery holder by sliding it out of the display assembly (C).
5. Replace the batteries with new ones. *
6. Carefully position the battery holder into the display assembly.
7. Connect the battery holder cable connector to the scale-battery connector and carefully slide the battery holder into the display assembly.
8. Attach the battery cover to the display assembly and install the screw.

* We recommend the use of EVEREADY Energizer® e2™ batteries.

AC ADAPTER USE

The ProPlus® display has been designed to operate on 9V DC and ships with an AC Adapter (APDT31). To operate the display using the AC adapter, insert the adapter into the jack located on the rear of the display. Use only a Health o meter® Professional adapter with the ProPlus® display.

QUICK START INSTRUCTIONS

WEIGHING A PATIENT

1. Make sure there is no object on the weighing platform.
2. Press the ON/OFF of I/O button to turn the display on.
3. Wait until -0.0 lb and $-ZERO$ on the left side of the display appear.
4. Ask the patient to step on the scale. The display will read $-WEIGHING$ until the weight calculation is complete and the weight is displayed. Depending on the movement of the patient on the scale, it may take several seconds for the scale to lock onto the weight.
5. If you wish to reweigh without asking the patient to step off and on the scale again, press the REWEIGH button.
6. Ask the patient to step off the scale.

NOTE: The scale will always default to the settings and units last used.

WARNING: If the scale will not be used for a prolonged period of time, remove the batteries to avoid a safety hazard.

QUICK START INSTRUCTIONS (CONTINUED)

TARE FUNCTION

When using the scale, the weight of an object such as a wheelchair or shoes accompanying the patient can be subtracted from the total weight of the patient alone. The Tare function automatically performs this subtraction, avoiding the need for manual calculations. ProPlus® scales allow tare weight to be entered automatically (PUSH BUTTON TARE) or manually using the keypad to enter the exact tare value (KEYPAD TARE).

AUTOMATIC (PUSH BUTTON) TARE

The user can set a tare value by pressing the TARE button (9) while there is a weight on the scale platform. The display will reset to zero and the word "TARE" will be displayed to indicate there is a tare value stored in memory.

NOTE: Due to the scale's sensitivity, we recommend using the REWEIGH function prior to setting the TARE weight, in order to eliminate any operator interference with the item to be tared.

MANUAL (KEYBOARD / KEYPAD) TARE

1. Make sure there is no object on the weighing platform.
2. Press the ON/OFF or I/O button to turn the display on.
3. Wait until "0.0 lb" and "ZERO" on the left side of the display appear.
4. Press the TARE button (9). The user will be prompted to enter the TARE value.
5. Use the keypad to enter the weight of the object desired to be tared (i.e. wheelchair) and press ENTER.
6. The value entered will be displayed as a negative value.
7. Place the patient and the tared object on the scale. The display will automatically deduct the entered tared weight from the gross weight of the patient and tared object.
8. The weight of the patient will appear on the scale.
9. The tared value is stored in memory until it has been changed, cleared or the display has been turned off.

Important: The TARE weight cannot exceed:
250lb (4021, 4025, 4026, 4261, 4650 & 4700)
200lb (4024)
150lb (4011)
125lb (4020)

REMOVING THE TARE

Once a tare value has been stored and the word "TARE" appears on the display, the tare value can be deleted from memory by pressing the TARE button. The scale will return to normal operation. You can also press the EXIT button.

CALCULATING BMI

1. Complete steps 1 to 4 on page 7 for "Weighing a Patient".
2. Press the BMI button (2).
3. The scale will prompt you to enter the patient's height. Use the keypad to enter the height.

Note: When weighing in lbs., height is entered in ¼" increments. For the fractional portion of the height, press 1 for ¼", press 2 for ½" and 3 for ¾". When weighing in kg the height is entered in 1cm increments.


4. Press ENTER.
5. Display will show the patient's BMI.
6. Press ENTER to return to normal weighing operation.

Note: The scale will not calculate a BMI for a patient 26.4 lb/12 kg or less. If a height is not entered within 30 seconds of pressing the BMI key, the scale will return to normal operation.


OPERATION INSTRUCTIONS

MENU

In the menu screen the user can set preferences and/or instruct the scale how to handle stored data. The menu can be navigated using the up and down keys (▲▼) or by entering the associated menu position number with the keypad. The menu has a “roll-over” way of working: when the user scrolls to the bottom of the menu and presses the down button, it will return to the top of the menu.


NOTE: The default mode of the scale is set so the menu option is locked out. If the menu button is pushed before activating the menu option, “MENULOCK” will show on the display. To temporarily reactivate the use of the MENU key, press and hold the MENU key for 5 seconds. During the menu access delay, the menu will display “MENULOCK”. After 5 seconds, the menu will show on the display and you may begin navigating through the menu as follows:



01 WEIGHT DISPLAY MODE (only applies to pound values, NOT metric values)

Before making changes to this scale setting, please refer to “Menu” instructions above to temporarily activate the Menu key.

The user can set the screen display value in either fraction of pound ($\frac{1}{4}$, $\frac{1}{2}$ or $\frac{3}{4}$ lb), in pounds and ounces (resolution of 4 ounces) or in decimals (resolution of 0.2 lb). When kg is selected as the units of use, these settings have no effect. The mode that is selected is used through all the screens.

OPERATION INSTRUCTIONS

02 TARE WEIGHT



Before making changes to this scale setting, please refer to Menu instructions on page 7 to temporarily activate the Menu key.

In addition to the two methods of entering tare values mentioned in the Quick Start Instructions, a tare value can be stored into memory by selecting option 02-Tare Weight of the User Menu .

The Tare Weight window will appear and prompt the user to enter the tare value using the numbered keypad and to press ENTER.


NOTE: The tare weight must be entered using the following increments: 4oz, 0.2lb, ¼ lb.

When the tared weight is not on the scale, the value entered will be display as a negative value. After the TARE has been entered, the scale goes back to normal operation. This TARE value is stored in memory until changed, cleared or the scale is turned off. If the tare weight is not entered in 30 seconds the scale will revert back to normal operation, and the Menu Lock will be activated.

The number that is to be changed will flash and will move from the left to the right after the appropriate number was entered or by using the ▲ key (left) and the ▼ key (right).

The word "TARE" appears on the display to indicate a weight has been tared. When the tared weight is removed from the scale the value will be displayed as a negative value.

OPERATION INSTRUCTIONS

Removing the Tare


Once a tare value has been stored and the word —TARE” appears on the display, the tare value can be deleted from memory by pressing the TARE button. The scale will return to normal operation.

03 DATA MANAGEMENT


Before making changes to this scale setting, please refer to the Menu instructions on page 9 to temporarily activate the Menu key

The scale managed patient data including patient ID, weight, height, tare, and BMI. The tare is stored in memory or transferred to a PC. This function is performed by opening a new data file. The maximum capacity of the scale is 270 files of different data.

1. Press the ID button (7) See page 4 for the illustration
2. Using the keypad, type in the identification number
3. Press ENTER


01 Auto-Download


Automatic Download is the default option and will transfer the value to the computer as soon as the patient steps off the scale or when the user presses the HOLD/RELEASE button if it was kept in —HOD”. If no computer is connected, the value is not transferred or stored and will be lost after the load is removed from the scale.


02 Store in memory

The value is stored in memory for later download to a computer. If the memory is full the user will be warned and given the option to transfer all values to the computer or to clear the memory of values. Press EXIT to return to the previous screen.


03 Do not store data

No data will be stored. Press EXIT to return to the previous screen.


OPERATION INSTRUCTIONS

04 Transfer now


All the values stored in memory are transferred to the computer and the scale memory is cleared of all values. If the transfer was unsuccessful, the values are kept in memory until successfully transfer or cleared. Press EXIT to return to the previous screen.

05 Clear memory


If Yes is selected and Enter is pressed, all the values in memory will be cleared. If No is selected, the display will return to the last menu. Press EXIT to return to the previous screen.

04 SCALE SETTINGS

Before making any changes to this scales settings, please refer to the Menu instructions on page 7 to temporarily activate the Menu key.


01 Auto Hold Time

The user can determine how long to display the weight reading once it is determined, regardless of whether the patient remains on the platform. The scale defaults to no Auto Hold Time. The maximum setting is 20 seconds Hold Time. To set the time, use the keypad to enter the desired number of seconds, and then press the ENTER key. The display will then return to the Scale Settings menu.

OPERATION INSTRUCTIONS


02 Sleep Time

The user can set the time elapsed before the scale goes into the sleep mode. The default is 1 minute. When the scale goes into sleep mode, STANDBY is displayed on the screen. The maximum setting is 10 minutes Sleep Time. Use the keypad to enter the desired amount of time and then press the ENTER key.


03 Auto Off Time


The user can determine how long the scale will operate before turning off automatically due to inactivity. Default time is 10 minutes. If the value is set to zero, the auto off function is disabled. The maximum setting is 60 minutes Auto Off Time. Use the keypad to enter the desired amount of time and then press the ENTER key.

04 Tone Volume


There is an option to adjust the beeping tone of the scale. This tone should sound when the scale has determined the weight on the platform, when a key is pressed, after the scale is turned on, at the end of self-test, or in the case of fault or warning. Press the ENTER key to set the volume. Press the EXIT key to return to the Scale Settings menu.


Use the ▲ and ▼ keys on the keypad to adjust the volume. Whenever the user presses the key to change a volume, a beep will sound to indicate the set volume level. If 0 is selected, the tone is disabled.

OPERATION INSTRUCTIONS

05 Disp Date & Time


06 Display Backlight

The user can set the brightness of the backlight. After selecting the desired brightness level, press ENTER to store the brightness level and to revert back to the Scale Settings menu.


07 Display Contrast

The user can set the contrast of the LCD. After selecting the desired contrast level, press ENTER to store the contrast level and to revert back to the Scale Settings menu.

OPERATION INSTRUCTIONS

08 Live Weight

By selecting “Yes” the user deactivates the motion-sensing mode and activates the Live Weight mode. In Live Weight mode the weight displayed will fluctuate with the patient’s movement. In Live Weight mode the scale will not lock on to the weight as it does when the motion-sensing mode is active.


To use motion sensing mode to verify a weight determined in the Live Weight mode, press the keypad REWEIGH button and the weight will be displayed on the screen.

To deactivate Live Weighing mode, highlight “No” and press enter.

09 Exit

This will take you out of the menu to the actual weight display.


05 SYSTEM SETTINGS

Before making any changes to this scales settings, please refer to the Menu instructions on page 7 to temporarily activate the Menu key.


01 Set Time & Date

The user can set the time and date using the keypad.


OPERATION INSTRUCTIONS

01 Set Time



To set the time move between hours, minutes and seconds using the up and down keys (▲▼) and enter the values on the keypad. After setting the time, press ENTER to select AM or PM. Press ENTER to save.

02 Set Date



Set the date using the up and down keys (▲▼) to move between month, day and year, and enter the values on the keypad. Press ENTER to save.

02 About


This screen displays the software version of the display. Either EXIT or ENTER will return to the previous menu.

03 Display 24 Hour

There are two options to display time in 12 hour (AM/PM) or 24 hour format


24-hour format

12-hour format


OPERATION INSTRUCTIONS

To display time in a 24 hour format, highlight “Yes” and press the ENTER button. The display will return to the System Settings menu. If the 24 hour format was selected and the correct time was set in the 12 hour format, the display will automatically convert the time to the 24 hour format. If the user needs to set the time and date, select menu option 01-Set Time and Date to set the correct time for the 24 hour format.


06 SYSTEM TEST

Before making any changes to this scales settings, please refer to the Menu instructions on page 7 to temporarily activate the Menu key.

01 Battery Test

When there are batteries in the scale, the display will show the estimated amount of battery life remaining until the batteries will have to be replaced.

NOTE: In order to complete the battery test, the display must be powered by batteries only. Unplug the display from the AC power


OPERATION INSTRUCTIONS

02 USB Connection


The display will test the connection to the computer and will show a message —Connection is OK” or —NO Connection”. If —NO Connection” is shown, check your USB connections to the display and your computer and retest the connection. Refer to qualified service personnel if problem persists.


03 Diagnostic Test

The display has a diagnostic routine where it tests the User Interface (UI) hardware functionality (LCD, keypad). In order to do this the user has to press all the keys according to the messages displayed on the screen.


If the requested command was not received or the wrong button was pressed, the following message will be displayed.

If after 10 seconds the requested command was not received, the following message will be displayed. If —UI Error Failure” is displayed, refer to qualified service personnel.

07 USER SETTINGS

Before making any changes to this display’s settings, please refer to Menu instructions to temporarily activate the Menu key.

OPERATION INSTRUCTIONS

01 Retain Entered Values

01-Retain Ent Values


Yes
No

▼ ▲ to change
Enter

This option allows the user to use the same values for ID, height and TARE between weighings. If this option is disabled, the user has to re-enter these values for each reading. If the values are not entered, only the weight is stored.

NOTE: These values cannot be retained by ID number.

02 Display Height & ID

When the user selects to display the height and ID of the patient, it will be displayed at the bottom of the screen. We recommend the use of this function to ensure that the patient's correct ID and height have been entered.

03 Prompt For Height

03-Prompt For Height

Yes
No

▼ ▲ to change
Enter

When this option is activated, the user will be asked to enter the patient's height after every weighing. The operator has 30 seconds to enter height.

10/14/2009
11:30pm

90.9 lb

Please enter Height:
0ft00 0/0in

▼ ▲ to change
Enter

04 Prompt For ID

04-Prompt For ID

Yes
No

▼ ▲ to change
Enter

When this option is activated, the user will be asked to enter the patient's ID after every weighing. The operator has 30 seconds to enter ID.

10/14/2009
11:30pm

90.0 lb

Please enter ID:
234587336

▼ ▲ to change
Enter

MAINTENANCE

GENERAL

This section provides instructions for maintenance, cleaning, troubleshooting and operator replaceable parts for ProPlus® Display Modules. Maintenance operations other than those described in this section should be performed by qualified service personnel.

MAINTENANCE

Before first use and after periods of non-use. Check display for proper operation and function. If the display does not operate correctly. Refer to qualified service personnel.

1. Check overall appearance of the display for any obvious damage, wear and tear.
2. Inspect the AC Adapter for cord cracking or fraying or for broken or bent prongs.
3. See parts list on page 22 if replacement parts are required.


CLEANING

Proper care and cleaning is essential to ensure a long life of accurate and effective operation.

Disconnect the scale from the AC power source.

1. Clean all external surfaces with a clean damp cloth or tissue. Anti-bacterial soap and water may be used. Dry with a clean soft cloth.
2. Do not immerse the display into cleaning or other liquid solution.
3. Do not use Isopropyl Alcohol or other cleaning solutions to clean the display.
4. Do not use abrasive cleaners.

EXPLODED DIAGRAM– DISPLAY MODULE

DISPLAY MODULE PARTS LIST

Field Replaceable parts				
Key no.	Part no.	Description	Qty.	REMARKS
	B1109401-0	Display Assembly	1	Includes all parts (Key no. 1-22).
3	B2033801-0	BATTERY HOLDER. PP	1	
12	B3822801-0	BATTERY COVER	1	
Parts listed below are for reference only (refer to display diagram on page 21)				
1		KEYPAD	1	
2		LENS, TRANSPARENT PC	1	
4		STATIC PROTECTION MEMBRANE	1	
5		LCD PCB ASM	1	
6		WIRE	1	
7		PCB ASM	1	
8		STATIC PROTECTION MEMBRANE	1	
9		SINGLE FACE GLUE FOAM	1	
10		INPUT/ OUTPUT POSTS LABEL	1	
11		SCREW M3X6	2	
13		SCREW M3.5X8	6	
14		HEAD BACK, ABS+Q235A	1	
15		SCREW, ϕ 3X6	4	
16		TERMINAL	1	
17		SCREW, ϕ 2.3 X 8	4	
18		WIRE	1	
19		POSITIVE PLATE	1	
20		NEGATIVE PLATE	1	
21		POSITIVE AND NEGATICE PLATE	1	
22		HEAD FRONT, ABS+Q235-A	1	

CALIBRATION PROCEDURE

All Health o meter[®] Professional scales are factory calibrated and do not require calibration prior to use.

Calibration can be performed in kg or lb, according to the units of measure used upon entering into the calibration mode.

Operator Action	Display
1. Make sure the display is turned off, then press and hold the HOLD/RELEASE button while pressing the ON/OFF or I/O button.	Enter load weight: 400.0LB (or 200.0KG)
2. Using the ▼ and ▲ keys and/or the numeric keypad set the calibration weight to at least 1/2 of the capacity of the scale (for the best results and greater accuracy, use the full capacity of the scale).	Enter load weight: 400.0LB (or 200.0KG)
3. Make sure there is nothing on the scale, then press ENTER button.	Zero calibration: Please clear the scale
4. Please wait 1-3 seconds without touching or shaking the scale for Zero Calibration to finish.	Zero calibration: Please wait
5. Load the scale with the required weight and press the ENTER button.	Weight calibration: Put: XXX.X
6. Please wait 1-3 seconds without touching or shaking the scale for Calibration to finish.	Weight calibration: Please wait
7. If calibration is successful, the display will show Factor OK. If the process failed, the display will show Factor Fail.	Factor Update Factor OK! (or Factor Fail!)
8. If successful, the display will begin to reboot.	Rebooting... Please wait
9. The display will return to the start up message.	Health o meter Pro-Plus
10. Wait until the display resumes normal operation and then remove the weight from the scale.	Health o meter Pro-Plus

PC COMMUNICATION PROTOCOL

This ProPlus[®] display uses an escape protocol for communicating with a PC and over the USB and serial ports. In an Escape protocol, the escape (<esc> or ASCII 27) is used to indicate that there is a command, and not data, following.

Two scenarios have been defined: Scale initiated communication and PC initiated communication. The following is a table of what can be sent across the communications line.

Scale Initiated	ESC Value	PC Initiated	ESC Value
Send Single reading	R	Update Firmware	U
Send continuous reading	R	Request current values/settings	Q
Send bulk readings	R	Diagnostics	A
Send Diagnostics	A	Send scale control messages	C
		Request stored data	R
		Set stored data	S
		Delete stored data	X

PC COMMUNICATION PROTOCOL

The following is a complete list of ESC characters that will be used:

Name	ESC character	ESC value with parameters	Description
Reading	R	R	This will tell the PC that the scale is sending a reading. Immediately following this will be the value that is sent (i.e. <ESC>R<ESC>I1234567890<ESC>W200.00<ESC>Nm<ESC>...<ESC>E
ID	I	Iccccccccc	This is the patient ID. (10 chars)
Weight	W	Wnnn.nn	This is the patient weight (i.e. W200.05 means 200.5)
Height	H	Hnnn.nn	This is the patient height
TARE	T	Tnnn.nn	This is the TARE weight that the scale is taring out
BMI	B	Bnnn.n	This is the patient's BMI
End Of Packet (EOP)	E	E	This indicates that the end of the command/datapacket has been reached
Units	N	Nc	This indicates in which unit the values have been taken (m=metric, c=constititional)
Power Status	J	J	Requests the power status of the scale. The response to this will be either a <ESC>O<ESC>E or <ESC>F<ESC>E or <ESC>L<ESC>E
On	O	O	When prompted for a power status, this will indicate that the scale is turned on
Off	F	F	When prompted for a power status, this will indicate that the scale is turned off (or is in low power/sleep mode) When in Low Power mode, the scale will return L
DateTime	D	Dnnnnnnnnnn n	This is a time and date string (MMDDYYYYhhmmss)
Update	U	Uc	This command has a parameter c, which will help doing the update of the firmware. First of all the PC will initiate the update by sending an Ur (Request for update). The scale will reply with Uc (Clear to send) or Ue (error, don't send). Then the scale will send U127K or U127B (the length of the data that is following in Kbytes or Bytes) and send the data.
Request	Q	Qccc	This is a request to receive scale global values such as brightness of LCD, auto hold time, keypad beep volume, etc. (See Table A) (i.e. <ESC>QAHT will return the Auto Hold Time of 10 seconds like this: <ESC>PAHT=10<ESC>E)
Response	P	Pccc=ccccccc	This will be the scales response to the request mentioned above
Diagnostics(1)	A	Accc	This is the request for a diagnostic test on certain parts of the scale (such as battery life, loadcells, keypad, LCD etc.). (See Table B)
Diagnostics(2)	Z	Zccc	This will be the response of the diagnostics done on the scale. Values will include error codes to indicate what is wrong with the scale, or all zeros to indicate that all is well.
Control (set a value)	C	Cccc=ccccccc	This is to set the value of the scale's global settings such as backlight, auto hold time, keypad beep volume etc. (See Table A) (i.e. <ESC>CAHT=20<ESC>E will set the Auto Hold Time to 20 seconds)
Get subject record data and readings	G	G	This will return all the information from both the subject record data as well as the readings stored in the readings table.
Set subject record data	S	S	This will set the subject record data for a specific patient. Immediately following this will be the value that is set. (i.e. <ESC>S<ESC>I1234<ESC>T200<ESC>Um<ESC>E will set the TARE for patient ID1234 to 200kg)
Clear readings	X	X	This will clear all the values stored in the readings table on the scale. Only the subject record data will not be deleted.
Clear subject record data	Y	Y	This will clear the subject record data for a specific patient. The patient ID will follow the Y immediately (<ESC>Iccccccccc). If no ID is given, ALL the data will be cleared.

We used descriptive characters to show the type of the parameters. These are: c=Character, n=Numeric. The number of parameters indicates the set number of placeholders for these characters and numbers (i.e. Zccc indicates that a value containing 3 characters will be given.)

PC COMMUNICATION PROTOCOL

Name of Control	Identifier	Unit
Auto Hold Time	AHT	nn (Seconds)
Auto Off Time	AOT	nn (Minutes)
Sleep Time	STM	nn (Seconds)
Keypad Volume	VKE	nn (01..10)
Stable weight Volume	VSW	nn (01..10)
Alarm Volume	VAL	nn (01..10)
Display of Date/Time	DAT	nn (01 or 10)
Backlight	BAC	nn (01..10)
Contrast	CON	nn (01..10)
Weight Display Mode (fractions(f), decimals(d) or Pounds and ounces (p))	WDM	c (f,d,p)
Tare Weight	TAR	nnn.nn (Kg)
Unit Of Measure (Metric or constitutional)	UOM	c (m or c)
Positions left in RAM	RAM	nnnn (0..max)

Table A - Scale Global Values List and Identifiers

Name of component to test	Identifier	Error Code
Load Cell A	LCA	E3A
Load Cell B	LCB	E3B
Load Cell C	LCC	E3C
Load Cell D	LCD	E3D
Battery	BAT	E4L (Bat Low, but still usable) or E4U (Bat Low and Unstable)
PC Communication (USB)	PCC	E05
Write to RAM	WRM	E08
Read from RAM	RRM	E09

Table B - Components to test

The different tables on the scale:

There will be two data tables on the scale. The first will store the subject record data (patient-related data that does not change) and the second will store the data recorded during the readings made.

Subject record data


Fields: *ID, TARE, HEIGHT*

Readings


Fields: *Index (referring to the ID from the Subject Record Table), Weight, Date-time*

Sample and explanation of esc protocol

1. When the scale measures the weight of a patient and sends this over the communications line to the PC, it will look like this:


PC COMMUNICATION PROTOCOL


2. When the PC requests the values of the stored records, it will look like this:

First the PC will ask the scale to send the data. We do this with the "Get Subject Record Data and Readings" command.

Then the scale will reply with the stored data.

If no data is stored on the scale, it will simply send


<ESC>R<ESC>E

PC COMMUNICATION PROTOCOL


Pin Configuration for ProPlus® Cables

The pins of the two connectors are connected as follows:

RJ45	function	DB9
7	TX	2
5	RX	3
8	GND	5

RJ45

DB9

Transmission Parameters

Baud rate	9600
Parity	None
Databits	8 bits
Stopbits	1bit
Startbits	1 bit
Hardware handshake	None

WARRANTY

3-YEAR LIMITED WARRANTY

What does the Warranty Cover?

Health o meter® Professional ProPlus® scales are warranted from date of purchase against defects of materials or in workmanship for a period of three (3) years. If you have an equipment problem that you cannot resolve, and the Customer Service determines that a defect of material or in workmanship exists, the customer's sole remedy will be repair or replacement of the scale at no charge. Replacement will be made with a new or remanufactured product or component. If the product is no longer available, replacement may be made with a similar product of equal or greater value. All parts including repaired and replaced parts are covered only for the original warranty period.

Who is Covered?

The original purchaser of the product must have proof of purchase to receive warranty service. Please save your invoice or receipt. Authorized dealers or retail stores selling Pelstar products do not have the right to alter, or modify or any way change the terms and conditions of this warranty.

What is Excluded?

Your warranty does not cover normal wear of parts or damage resulting from any of the following: negligent use or misuse of the product, use on improper voltage or current, use contrary to the operating instructions, abuse including tampering, damage in transit, or unauthorized repair or alternations. Further, the warranty does not cover Acts of God, such as fire, flood, hurricanes and tornadoes. This warranty gives you specific legal rights, and you may also have other rights that vary from country to country, state to state, province to province or jurisdiction to jurisdiction.

To get Warranty Service make sure you keep your sales receipt or document showing proof of purchase. Call the Customer Service Center nearest you for assistance-- please see page 6 of this manual. **Before returning a product, you must obtain authorization from Customer Service.** Customer Service personnel will give you a return authorization number. Please note this number on the outside of your shipping box. Returns without a Return Authorization number will not be accepted for delivery. **Attach your proof of purchase to your defective product along with your name, address, daytime telephone number and description of the problem. Carefully package the product and send it with shipping and insurance prepaid to the address you are given by Customer Service.** Write, "Attention RA# _____", and the shipping address, on the shipping label, following Customer Service instructions.

PRODUCT LIFE: 10 YEARS / 100,000 CYCLES

ProPlus® is a registered trademark of Pelstar, LLC. Health o meter® is a registered trademark of Sunbeam Products, Inc. used under license. Health o meter® Professional products are manufactured, designed and owned by Pelstar LLC

We reserve the right to improve, enhance or modify Health o meter® Professional product features or specifications without notice.

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© Pelstar, LLC 2011