# Reichert® UNISTAT® Bilirubinometer

Model 1310310C/1310311C

**Instruction Manual**

**Table of Contents**

<b>Section</b>	<b>Page</b>
1.0 APPLICATION .....	4
2.0 INSTALLATION .....	4
3.0 PRINCIPLES OF OPERATION .....	4
4.0 PERFORMANCE CHARACTERISTICS AND SPECIFICATIONS .....	6
4.1 Sample Type .....	6
4.2 Specifications .....	6
5.0 OPERATING INSTRUCTIONS .....	8
5.1 Ready Mode .....	8
5.2 Preparation for Use .....	8
6.0 CALIBRATION .....	12
6.1 Calibration Frequency .....	12
6.2 Calibration Material .....	12
6.3 Calibration Procedure .....	12
6.4 Calibration Verification .....	14
7.0 QUALITY CONTROL .....	14
7.1 Use of Commercial Serum Controls .....	14
7.2 Quality Control Frequency .....	14
8.0 OPERATIONAL PRECAUTIONS AND LIMITATIONS .....	16
8.1 Specimen Collection and Preparation .....	16
8.2 Known Interfering Substances .....	16
8.3 Recommended Storage & Handling of Specimens .....	16
9.0 RESULTS .....	16
10.0 EXPECTED VALUES .....	16
11.0 SERVICE .....	18
11.1 Bulb Replacement .....	18
11.2 Adjusting Calibration Value .....	18
11.3 Changing Units (mg/dL — µmol/L) .....	20
11.4 Fuse .....	20
11.5 Care of the Instrument .....	20
12.0 HAZARDS .....	20
13.0 ACCESSORIES .....	22
14.0 BIBLIOGRAPHY .....	22

**Explanation of Warning Symbols**

**CAUTION:** Replace with same type and rating fuse

1310310C: T 600ma, 250V, 5x20mm  
1310311C: T 250ma, 250V, 5x20mm

**CAUTION:** Risk of electric shock

**CAUTION:** (Refer to accompanying documents)

**NOTE:** This manual is for use with the 1310310C/1310311C REichert UNISTAT Bilirubinometer.

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 **ATTENTION:** (Se reporter aux documents inclus)

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## 1.0 APPLICATION

The REICHERT UNISTAT® Bilirubinometer Model 1310310C/1310311C is a direct spectrophotometer for measuring total bilirubin in neonatal serum to aid in indicating the risk of bilirubin encephalopathy (kernicterus). The REICHERT UNISTAT Bilirubinometer is intended for invitro diagnostic use only.

**CAUTION: THIS INSTRUCTION MANUAL AND THE FIGURES CONTAINED HEREIN DESCRIBE THE REICHERT UNISTAT BILIRUBINOMETER (MODEL 1310310C, 1310311C) AS CLOSELY AS POSSIBLE. PLEASE CONTACT US AT (716) 686-4500 IF YOU HAVE DIFFICULTY UNDERSTANDING THE MANUAL OR INSTRUMENT USE.**

**RESPONSIBILITY: IT IS THE RESPONSIBILITY OF THE USER TO OPERATE THE INSTRUMENT PROPERLY ACCORDING TO THE INSTRUCTIONS IN THIS INSTRUCTION MANUAL; INCORRECT OPERATION MAY INVALIDATE RESULTS. DIAGNOSES AND CONCLUSIONS DERIVED FROM DATA OBTAINED WITH THE INSTRUMENT ARE NOT THE RESPONSIBILITY OF THE MANUFACTURER.**

**CAUTION: NEVER USE A TWO PRONG ADAPTER BETWEEN THE POWER CORD AND THE POWER SOURCE-ELECTRICAL SHOCK HAZARD MAY RESULT.**

## 2.0 INSTALLATION

The REICHERT UNISTAT Bilirubinometer is shipped ready for installation. Be sure to check the data plate on the rear of the unit for the appropriate power source. Install instrument on a hard level surface.


The power switch is located on the back of the instrument, near the data plate.

When first powered up, the instrument will briefly display its revision number and "LLL" indicating that a system check is in progress.

After a five-minute warm-up period, the instrument is ready for calibration (see Section 6.0). Calibrate the instrument, verify calibration and run quality control samples before analyzing patient samples.

## 3.0 PRINCIPLES OF OPERATION

The REICHERT UNISTAT Bilirubinometer is a direct spectrophotometer that measures total bilirubin in neonatal serum samples; dilution is not necessary. Light is passed through the sample cuvette. The beam is split by a dichroic mirror and directed toward the 460nm and 550nm filters. The transmittance of the light exiting the filters is measured by the photodetectors. The control electronics calculate

**Figure 1** Optical Schematic for Reichert UNISTAT Bilirubinometer  
**Abb. 1** Schaubild der Optik des Reichert UNISTAT Bilirubinometer  
**Figure 1** Schéma optique des Reichert bilirubinomètres UNISTAT  
**Figura 1** Esquema óptico del bilirubinómetro UNISTAT Reichert

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total bilirubin based on output from the photodetectors (see Fig. 1).

Since peak absorbance is shared by bilirubin and oxyhemoglobin at 460nm, it is necessary to measure oxyhemoglobin at 550nm (see Fig. 2). The difference in the readings at 460nm and 550nm is the total bilirubin value, corrected for the presence of oxyhemoglobin.

Since the Reichert UNISTAT Bilirubinometer measures and subtracts oxyhemoglobin, bilirubin results are typically not affected by hemolyzed samples.

#### 4.0 PERFORMANCE CHARACTERISTICS AND SPECIFICATIONS


##### 4.1 Sample Type

The REICHERT UNISTAT Bilirubinometer is designed to measure total **bilirubin only in patients who are classified as neonates by the clinician**. Interfering substances become more prevalent in patients who are not neonates (see Section 8.0).

##### 4.2 Specifications

- Sample size:** 20µL
- Measuring range:** 0-40 mg/dL (0-684 µmol/L)
- Results displayed:** 5 sec.
- Turnaround time:** 15 sec.
- Accuracy:** ±5% (0-30 mg/dL/0-513 µmol/L)  
±10% (31-40 mg/dL/530-684 µmol/L)
- Precision:** Precision was evaluated as outlined in NCCLS Protocol EP5-T2 (see Sec. 14.0, Ref. 1).

Level	Mean (mg/dL/µmol/L)	Precision Within Run (% C.V.)	Total (% C.V.)
I	1.2/21	10	15
II	7.6/130	5	10
III	22.2/380	5	10

**Figure 2** Absorption of Bilirubin and Oxyhemoglobin  
**Abb. 2** Lichtabsorptionsgrad von Bilirubin und Oxyhämoglobin  
**Figure 2** Absorption de la bilirubine et de l'oxyhémoglobine  
**Figura. 2** Absorción de bilirrubina y oxihemoglobina

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## Comparison of Methods:

The REICHERT UNISTAT Bilirubinometer and the Jendrassik-Grof (Dumas) proposed reference method were compared for total neonatal bilirubin. A combination of clinical samples from normal and hyper-bilirubinemic neonates and spiked serum control samples were analyzed following NCCLS Protocol EP9-P (Figure 2b) (see Section 14.0, Ref. 2).

**Results displayed: 5 sec.**

Regulatory compliance: The Reichert UNISTAT Bilirubinometer is registered and listed with Intertek and carries the following marks: cETLus, CE.

- Model 1310310C: 115V AC, 60Hz, 60W
- Model 1310311C: 230V AC, 50 Hz, 60W
- Indoor use only
- RH 80% up to 31°C decreasing linearly to 50% RH at 40°C
- Pollution degree 2
- Altitude up to 2,000m
- Temperature 5°C to 40°C
- Mains supply voltage  $\pm 10\%$
- Installation CAT II

## 5.0 OPERATING INSTRUCTIONS

### 5.1 Ready Mode

The REICHERT UNISTAT Bilirubinometer may be left in the standby (on) mode; component life will not be shortened by keeping the instrument in this mode. The standby mode is indicated by a lighted dot in the upper left corner of the LED display.

**Once an hour the bilirubinometer will perform a system check. During the system check, “LLL” will be displayed briefly.**

There is an on-off switch on the rear panel near the data plate which can be used to turn the power on and off.


### 5.2 Preparation for Use

Refer to Section 8.0 for sample handling considerations.

Before analyzing patient samples, calibration verification (see Section 6.4) and quality control procedures (see Section 7.0) **must be completed**. If calibration verification indicates the instrument is out of calibration, follow the procedure in Section 6.0 to recalibrate.

1. Using a pipette or syringe, **completely fill the cuvette** (Reichert cat. no. 13102140) with serum, holding the narrow portion of the cuvette between your thumb and index finger. The cuvette holds approximately 20  $\mu\text{L}$  of serum. Since the exact volume is determined by the cuvette, it is not necessary to measure the sample.

**IMPORTANT: ONLY REICHERT CATALOG NUMBER 13102140 DISPOSABLE CUVETTES MAY BE USED WITH THE UNISTAT BILIRUBINOMETER.**

**Figure 2b** Comparison of Methods  
**Abb. 2b** Die Methoden im Vergleich  
**Figure 2b** Comparaison des méthodes  
**Figura 2b** Comparación de los métodos

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**13102140 CUVETTES ARE FOR ONE-TIME USE ONLY AND NEVER SHOULD BE REUSED. USE OF OTHER CUVETTES, OR REUSED CUVETTES, MAY AFFECT REICHERT UNISTAT RESULTS.**

To avoid erroneous readings, the outside surface of the filled cuvette must be free of serum, dust, dirt and smudges; clean as necessary with a disposable wipe or gauze.

The sample seen in the clear rectangular area of the cuvette must be free of bubbles. If bubbles are observed, tap the cuvette lightly to clear them and refill as necessary.

2. Open the black door on the top of the instrument and insert the cuvette into the sample carrier until the widest portion of the cuvette is seated. Close the door (Fig. 3).

3. Press the START button. There will be a 2 second delay to allow the lamp to come to equilibrium, after which the sample carrier will rotate to bring the sample cuvette into the optical path. After a short pause, bilirubin concentration will be displayed.


4. After 15 seconds, the measurement cycle will be complete and the sample carrier will return the sample to the load position. Remove the cuvette. The display will show the bilirubin value for approximately 10 seconds, then clear.

If the measured value of the sample exceeds 40 mg/dL (684  $\mu\text{mol/L}$ ), the display will flash the result to warn of the change in the result's accuracy (see Section 4.0). If the result is 51.9 mg/dL (887  $\mu\text{mol/L}$ ) or higher "----" will appear on the display.

If the REICHERT UNISTAT Bilirubinometer is unable to measure the sample, one of the following error codes will appear:

Error Code	Problem	Corrective Action
E10	Cuvette not in sample path	<ul style="list-style-type: none"> <li>Turn instrument off and on, cuvette will rotate back to initial position. Seat cuvette completely in holder.</li> <li>Motor problem; Contact service.</li> </ul>
E11	Insufficient light	<ul style="list-style-type: none"> <li>Lamp is out; replace lamp.</li> <li>Light path blocked; contact service.</li> <li>Optical misalignment; contact service.</li> </ul>

5. After the display clears the previous value, the next sample may be run. The REICHERT UNISTAT Bilirubinometer does not require any rinsing or other maintenance between samples except for normal quality control and calibration procedures (see Sections 6.0 and 7.0).

**Figure. 3** Insert the cuvette into the sample carrier

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## 6.0 CALIBRATION

### 6.1 Calibration Frequency

The REICHERT UNISTAT Bilirubinometer requires calibration upon installation and then at least at six-month intervals, unless quality control data indicate the need for earlier calibration. The REICHERT UNISTAT Bilirubinometer must be recalibrated after all service procedures, including bulb changes.

### 6.2 Calibration Material


The REICHERT UNISTAT Bilirubinometer is calibrated with a stable assayed glass calibration cuvette (Reichert cat. no. 13103410) which is included with the instrument (Fig. 4). No other calibration material is to be used for calibrating the REICHERT UNISTAT Bilirubinometer. Replacement calibration cuvettes (Reichert cat. no. 13103410) may be ordered from your Reichert representative.

### 6.3 Calibration Procedure

1. Remove any smudges or dirt from the calibration cuvette using a soft material such as gauze. If necessary, clean the glass surface of the cuvette with a cotton swab dipped in methanol. Hold the cuvette by the narrow portion to prevent recontamination.
2. To begin the calibration cycle, depress and hold the START button until "ccc" appears on the LED display. Release the START button.
3. The display now will flash the value of the expected calibration cuvette (the instrument is set with the value of its calibration cuvette at the factory). The top of the calibration cuvette is labeled with its value, in mg/dL and  $\mu\text{mol/L}$  units.

**CAUTION: IF THE VALUE ON THE CALIBRATION CUVETTE DIFFERS FROM THE VALUE BEING DISPLAYED BY THE REICHERT UNISTAT BILIRUBINOMETER, DISCONTINUE THE CALIBRATION PROCEDURE BY TURNING THE POWER OFF AND ON (TO CLEAR) AND REFER TO SECTION 11.2.**

4. Open the black door and insert the specified calibration cuvette into the sample carrier until the widest portion of the cuvette is seated. Close the black door.
5. Depress and release the START button. The instrument will begin its calibration cycle. When the calibration cuvette returns to its original position, the calibration procedure is complete. E 12 will appear if the calibration cuvette measurement is more

**Figure 4** Calibration Cuvette and High-Level Check Cuvette

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than 1 mg/dL (17.1  $\mu\text{mol/L}$ ) different than the expected value. To correct:

- Be sure the correct cuvette (13103410) was used.
- Be sure the calibration value programmed matches the 13103410 calibration cuvette used (see Section 6.3).
- If the correct cuvette and value were used, the error may be the result of an optical problem. Contact your Reichert service representative.

E13 will appear if there is an electronic failure.

6. Remove the calibration cuvette.

7. See Section 6.4 for the calibration verification information.

A summary of the calibration error codes is shown in the following table:

Error Code	Problem	Corrective Action
E12	Standard reading incorrect	<ul style="list-style-type: none"> <li>• Clean standard</li> <li>• Verify (13103410) Cuvette and programmed calibration value</li> <li>• Optical problem: contact service</li> </ul>
E13	Calibration Failure	<ul style="list-style-type: none"> <li>• Board Failure: contact service</li> </ul>

## 6.4 Calibration Verification

Two assayed glass cuvettes are provided with the REICHERT UNISTAT Bilirubinometer:

1. Reichert cat. no. 13103410 Calibration Cuvette (assay values range from 19 to 23.9 mg/dL (325 to 409  $\mu\text{mol/L}$ ).
2. Reichert cat. no. 13103420 High-Level Check Cuvette (assay value approximately 40 mg/dL (684  $\mu\text{mol/L}$ )).

The cuvettes are labelled with actual assayed values.

Cat. no. 13103410 and 13103420 glass cuvettes may be used to check the bilirubinometer at the mid and high points of its 0 - 40 mg/dL (0 - 684  $\mu\text{mol/L}$ ) measuring range. A sample cuvette (Reichert cat. no. 13102140) filled with distilled water may be used to check zero.

The 13103430 Low-Level Check Cuvette (assay value approximately 12 mg/dL (205  $\mu\text{mol/L}$ )) may be purchased separately.

## 7.0 QUALITY CONTROL

### 7.1 Use of Commercial Serum Controls

Analysis of at least a normal and abnormal level of a commercial serum control, assayed for total bilirubin, is recommended for checking performance of the Reichert Unistat Bilirubinometer.

### 7.2 Quality Control Frequency

Quality control procedures should be performed before and after each sample run in response to questionable patient results and/or as mandated by local regulations.

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## 8.0 OPERATIONAL PRECAUTIONS AND LIMITATIONS

**IMPORTANT NOTE:** It is necessary to measure parameters other than total serum bilirubin concentration for assessment of the likelihood of kernicterus development, or clinical management of jaundiced neonates. A separate determination of "conjugated" bilirubin by conventional assay methods is also required.

### 8.1 Specimen Collection and Preparation

Neonatal serum samples may be collected using standard laboratory practice. Serum is added directly to the measuring cuvette; dilution is not required.

### 8.2 Known Interfering Substances

Pigments that absorb at 460nm, including, but not limited to lipochromes and carotenoids, will affect results. Therefore, only neonatal samples collected before these pigments are present may be accurately measured with the REICHERT UNISTAT Bilirubinometer.

Certain drugs, including, but not limited to amines, fatty acids, phenols and sulfonamides, may interfere with REICHERT UNISTAT Bilirubinometer measurements (see Section 14.0, Ref. 3).

Intravenous fat emulsions (including Lyposin®) have interfered with direct spectrophotometric measurements (see Section 14.0, Ref. 4).

Lipemic samples will interfere with REICHERT UNISTAT Bilirubinometer measurements.

### 8.3 Recommended Storage and Handling of Specimens

Under standard laboratory conditions, bilirubin samples are unstable to a degree varying from patient-to-patient. Samples must be protected from light to avoid photolysis of the serum bilirubin and be analyzed promptly (see Section 14.0, Ref. 5).

Serum stability is generally accepted as 1 week at 2°C – 8°C and 3 months frozen (-20° – 0°C) (R. Henry, "Clinical Chemical Principles and Technics", Harper & Row, Hagerstown, MD, 1974).

## 9.0 RESULTS

The REICHERT UNISTAT Bilirubinometer displays total bilirubin results directly from undiluted serum; calculations are not required.

### 10.0 EXPECTED VALUES

Normal serum neonatal bilirubin values are listed below. These values are provided only as a guide, rather than as absolute normals, and should be interpreted with discretion.

Age	Premature		Full-Term Newborn	
	mg/dL	µmol/L	mg/dL	µmol/L
Up to 24 hours:	1-8	17-137	2-6	34-103
Up to 48 hours:	6-12	103-205	6-10	103-171
3-5 days:	10-14	171-239	4-8	68-137

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(ver Sec. 14.0, ref. 6)

## 11.0 SERVICE

Contact your Reichert service representative for factory-authorized service information.

### 11.1 Bulb Replacement

The calculated lamp life for the REICHERT UNISTAT Bilirubinometer is over 10,000 measuring cycles.

**CAUTION: UNPLUG THE INSTRUMENT BEFORE PROCEEDING WITH BULB REPLACEMENT.**

When it is necessary to change the lamp (Reichert cat. no. 13110630), access the lamp through the panel (Fig. 5) in the bottom of the instrument by removing the two screws holding it in place. The lamp is a prefocused, double contact lamp bulb (GE 1460X). It is easily removed using the lamp extractor stored in the cover of the lamp compartment (see Fig. 6).

**RECALIBRATION OF THE REICHERT UNISTAT BILIRUBINOMETER WITH THE CAT. NO. 10341 CALIBRATION CUUVETTE, FOLLOWED BY CALIBRATION VERIFICATION AND QUALITY CONTROL PROCEDURES, MUST BE COMPLETED AFTER RELAMPING (see Sections 6.0 and 7.0).**

### 11.2 Adjusting Calibration Value

If the REICHERT UNISTAT Bilirubinometer displays a calibration value different from the value assayed on the calibration cuvette (during the calibration cycle), the instrument's expected calibration value must be adjusted. This procedure is typically required when a replacement calibration cuvette is obtained (the instrument is set with the value of its calibration cuvette at the factory).

**CAUTION: UNPLUG THE INSTRUMENT BEFORE PROCEEDING WITH SWITCH ADJUSTMENT.**

1. Access the switch bank through the panel (Fig. 5) in the bottom of the instrument by removing the two screws holding it in place.
2. Turn on switch 1 on the switch bank.
3. Plug the unit in and turn the power on.
4. Depress the start button (on the top of the instrument) until the value which matches the calibration cuvette's assayed value is displayed.
5. Return switch 1 on the switch bank to the "off" position.
6. Turn off and unplug the instrument. Reinstall the panel.
7. Power on instrument and proceed with calibration (see Section 6.0).

**Figure 5** Access Panel to Lamp and Switches

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**Figure 7** Switch Bank - Normal switch settings: units of measure in mg/dL: switches 1-6 off, switch 7 on units of measure in  $\mu\text{mol/L}$ : switches 1-7 off.

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### 11.3 Changing Units (mg/dL –µl/L)

**CAUTION: UNPLUG THE INSTRUMENT BEFORE PROCEEDING WITH CHANGING UNITS.**

1. Access the switch bank through the panel (Fig. 5) in the bottom of the instrument by removing the two screws holding it in place.
2. Switch 7 sets the unit of measure of the display. If this switch is in the "on" position, results will be displayed in mg/dL. If the switch is in the "off" position, units will be displayed in µmol/L.
3. If units of measure are changed from the factory setting, the unit label next to the LED display must be changed to correspond to the new units. Contact your Reichert service representative for a new overlay with the appropriate units.

**AFTER UNITS OF MEASURE ARE CHANGED THE CALIBRATION VALUE MUST *BE* RE-ADJUSTED** (see Section 11.2)

#### 11.4 Fuse

The REICHERT UNISTAT Bilirubinometer has two fuses located on the rear panel (Fig. 6). Use only the fuse specified on the data plate for replacement.

**CAUTION: UNPLUG THE INSTRUMENT BEFORE CHANGING ANY FUSES.**

#### 11.5 Care of the Instrument

The REICHERT UNISTAT Bilirubinometer can be cleaned with Naptha or soap and water. Avoid the use of ORA.

### 12.0 HAZARDS

**THE REICHERT UNISTAT BILIRUBINOMETER IS A LINE-OPERATED INSTRUMENT AND SHOULD BE CONNECTED TO A PROPERLY GROUNDED LINE. IT ALWAYS SHOULD BE DISCONNECTED FROM THE LINE BEFORE THE OUTER CASE IS REMOVED, AND BEFORE THE DOOR ON THE BOTTOM OF THE INSTRUMENT IS OPENED (FIG. 5).**

#### POTENTIAL BIOHAZARD

Identifies an activity or area where the operator may be exposed to potentially infectious materials or substances.

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### 13.0 ACCESSORIES

Cat. no.	Description
13102140	Disposable cuvettes (Box of 51)
13110630	Replacement lamp (Box of 3)
13103410	Calibration cuvette
13103420	High-level check cuvette
13103430	Low-level check cuvette

### 14.0 BIBLIOGRAPHY

1. Evaluation of precision performance of clinical chemistry devices. National Committee for Clinical Laboratory Standards (1992).
2. User comparison of quantitative clinical laboratory methods using patient samples (EP9-P). National Committee for Clinical Laboratory Standards (1986).
3. Young, D. S., Pesianer, L. C. and Gibberman, V., Effects of drugs on clinical chemistry tests. *Clin. Chem.* 21, special edition (April 1975).
4. Moore, John J., Sax, Sylvan M. and DeFranc, Susan, Lyposin® Interference with neonatal bilirubin measurements. *Clin. Chem.* 28, 2334-35 (1982).
5. Doumas, B. T. and Lemberger, M., The effect of temperature and light on the short-term stability of conjugated and unconjugated bilirubin solutions. *Clin. Chem.* 21 (Abstract), 1000 (1975).
6. Tietz, N. Fundamentals of Clinical Chemistry, Philadelphia: W.B. Saunders Co., pg 740 (1987).

#### The REICHERT UNISTAT Bilirubinometer is manufactured by:

Reichert Inc.  
3362 Walden Avenue  
Depew, NY USA 14043

Liposyn is a registered trademark of Abbott Laboratories, North Chicago, IL, USA.

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### 14.0 LITERATURVERZEICHNIS

1. Evaluation of precision performance of clinical chemistry devices. National Committee for Clinical Laboratory Standards (1992).
2. User comparison of quantitative clinical laboratory methods using patient samples (EP9-P). National Committee for Clinical Laboratory Standards (1986).
3. Young, D. S., Pesianer, L. C. and Gibberman, V., Effects of drugs on clinical chemistry tests. *Clin. Chem.* 21, Sonderausgabe (April 1975).
4. Moore, John J., Sax, Sylvan M. und DeFranc, Susan, Lyposin® Interference with neonatal bilirubin measurements. *Clin. Chem.* 28, 2334-35 (1982).
5. Doumas, B. T. and Lemberger, M., The effect of temperature and light on the short-term stability of conjugated and unconjugated bilirubin solutions. *Clin. Chem.* 21 (Abstrakt), 1000 (1975).
6. Tietz, N. Fundamentals of Clinical Chemistry, Philadelphia: W.B. Saunders Co., 740 (1987).

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### 14.0 BIBLIOGRAPHIE

1. Evaluation of precision performance of clinical chemistry devices. National Committee for Clinical Laboratory Standards (1992).
2. User comparison of quantitative clinical laboratory methods using patient samples (EP9-P). National Committee for Clinical Laboratory Standards (1986).
3. Young, D. S., Pesianer, L. C. and Gibberman, V., Effects of drugs on clinical chemistry tests. *Clin. Chem.* 21, édition spéciale (Avril 1975).
4. Moore, John J., Sax, Sylvan M. and DeFranc, Susan, Lyposin® Interference with neonatal bilirubin measurements. *Clin. Chem.* 28, 2334-35 (1982).
5. Dumas, B. T. and Lemberger, M., The effect of temperature and light on the short-term stability of conjugated and unconjugated bilirubin solutions. *Clin. Chem.* 21 (résumé), 1000 (1975).
6. Tietz, N. Fundamentals of Clinical Chemistry, Philadelphia: W.B. Saunders Co., 740 (1987).

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