

Smartphone Eye Imaging Adapter HSA-001

HSA001



Hans Heiss



OFFICIAL DISTRIBUTOR
MERCOFRAMES OPTICAL CORP

5555 NW 74 Ave. Miami, FL 33166. United States

Phone: 305-882-0120

sales@mercoframes.net

www.mercoframes.com

Product Introduction

QuikVue[®] filter kit includes two filter covers to work with QuikVue smartphone eye imaging adaptor. MeiboFilter[™] is a red filter for MGD screening. FluoFilter[™] is a yellow filter for fluorescein imaging.

VFR-100



Product Name: MeiboFilter[™]
Color: Red
Material: High transmittance glass
Application: Meibomian gland dysfunction screening
Imaging Mode: Phone camera in black/white mode
Image Processing App: Snapseed[®]

VFY-100



Product Name: FluoFilter[™]
Color: Yellow
Material: High transmittance glass
Application: Fluorescein imaging
Imaging Mode: Phone camera in color mode

*MeiboFilter[™] has only been tested on iPhone models so far. Its performance on other phone models hasn't been validated yet.

*Snapseed[®] is a registered trademark of Google.

202009 *Subject to change in design or specification without advance notification.

Hans Heiss



OFFICIAL DISTRIBUTOR
MERCOFRAMES OPTICAL CORP

5555 NW 74 Ave. Miami, FL 33166. United States

Phone: 305-882-0120

sales@mercoframes.net

www.mercoframes.com