

# Smartphone Eye Imaging Adapter

HSA001



Hans Heiss



**OFFICIAL DISTRIBUTOR**  
**MERCOFRAMES OPTICAL CORP**

5555 NW 74 Ave. Miami, FL 33166. United States

Phone: 305-882-0120

[sales@mercoframes.net](mailto:sales@mercoframes.net)

[www.mercoframes.com](http://www.mercoframes.com)

# Product Introduction

QuikVue<sup>®</sup> filter kit includes two filter covers to work with QuikVue smartphone eye imaging adaptor. MeiboFilter<sup>™</sup> is a red filter for MGD screening. FluoFilter<sup>™</sup> is a yellow filter for fluorescein imaging.

## VFR-100



Product Name: MeiboFilter<sup>™</sup>  
Color: Red  
Material: High transmittance glass  
Application: Meibomian gland dysfunction screening  
Imaging Mode: Phone camera in black/white mode  
Image Processing App: Snapseed<sup>®</sup>

## VFY-100



Product Name: FluoFilter<sup>™</sup>  
Color: Yellow  
Material: High transmittance glass  
Application: Fluorescein imaging  
Imaging Mode: Phone camera in color mode

\*MeiboFilter<sup>™</sup> has only been tested on iPhone models so far. Its performance on other phone models hasn't been validated yet.

\*Snapseed<sup>®</sup> is a registered trademark of Google.

202009 \*Subject to change in design or specification without advance notification.

**Hans Heiss**



**OFFICIAL DISTRIBUTOR**  
**MERCOFRAMES OPTICAL CORP**

5555 NW 74 Ave. Miami, FL 33166. United States

Phone: 305-882-0120

sales@mercoframes.net

www.mercoframes.com