



Sonomed Escalon™

Other Surgical Instruments

For Vitreoretinal Surgery

Iris Suture Device [TR9921]

Convenient tool for performing suture technique for iris retraction.

10 Devices per box, Sterile, Single-Use

Irrigating Contact Lens [TR9840-11; TR 9841-11]

Ideal for use during silicone oil procedures, irrigating contact lens enables visualization of surgical field while ensuring proper irrigation of cornea and orbit. Available in flat or minus styles, single-use lens is convenient requiring no oil clean up.

Packaged Individually, Non-Sterile (EtO-Compatible)

TR9840-11 Flat Lens

TR9841-11 Minus Lens

Retractable MVR Blade [TR1-0100]

Retractable MVR blade helps ensure needlestick safety in the OR and OSHA compliance. Blade can be locked in place in either engaged position for use during surgery or retracted position for pre- and post-use safety. Patented design allows single-hand operation to both engage and retract blade.

20 ga, Individually Packaged, Sterile, Single-Use



-
- 01 Iris Suture Device
 - 02 Irrigating Contact Lens
 - 03 Retractable MVR Blade
-



MERCOFRAMES OPTICAL CORP

5555 Nw 74 Ave. Miami, Fl. 33166. 305-882-0120

www.mercoframes.net

ale@mercoframes.net