

ADPTIK

A close-up photograph of an optical instrument, likely a lensometer or similar device. The image shows a lens being held in place by a metal frame. A scale is visible, with markings for 0, 10, 20, 30, 40, 50, and 60. The lens is partially visible, showing its curved surface. The background is a plain, light color.

TRIAL LENSES

Frames and lens sets



trial lens sets

precision set

Intuitive color-coded lens rings.

Paired Spherical Lenses

- ± 0.12
- ± 0.25 to $\pm 6.00D$ in $0.25D$ steps
- ± 6.50 to $\pm 8.00D$ in $0.50D$ steps
- ± 9.00 to ± 16.00 in $1.00D$ steps
- ± 18.00 , ± 20.00

Paired Cylinder Lenses

- ± 0.12
- ± 0.25 to $\pm 4.00D$ in $0.25D$ steps
- ± 4.50 ; ± 5.00 ; ± 6.00

Prism Lenses

- Paired 0.50 , 1 , 2 , 3 , 4Δ lenses
- Single 5 , 6 , 8 , 10Δ lenses

Auxiliary Lenses

- Red Lens
- Green Lens
- Slit Occluder
- $0.5mm$ Pinhole Occluder
- $1mm$ Pinhole Occluder
- Black Occluder
- White Maddox Rod
- Plano
- Cross Target




trial lens sets
precision set

Spherical Lenses

- ± 0.25 to ± 3.00 D in 0.25D steps
- ± 3.50 to ± 7.00 D in 0.50D steps
- ± 8.00 to ± 14.00 in 1.00D steps
- 16.00D

Cylinder Lenses


- 0.25 to 1.00D in 0.25D steps
- 1.50 to 5.00D in 0.50D steps



trial lens sets perimetry set

- Perimeter full aperture trial lens set with wooden case
- Available in plus or minus cylinder





trial frames

universal

- Make minute adjustments easily with the S4OPTIK universal trial frame
- Lightweight metal and plastic construction
 - Clear and precise scales



trial frames

fixed PD

- Achieve any correction with the essential fixed PD trial frame
- Each eye has 4 lens cells (3 in front and 1 behind)
- Available in 50, 52, 54, 56, 58, 60, 62, 64, 66 an 68mm sizes



Hand held cross cylinder

- Choose from hand held cross cylinders in ± 0.25 , ± 0.50 , ± 0.75 , ± 1.00 and ± 2.00 powers



trial clips

Halberg trial clips

- Featuring three cells for 38mm trial lenses and accessory lenses
- Helpful everyday and in low vision testing



MERCOFRAMES OPTICAL CORP.

📍 5555 NW 74 AVE. Miami, FL 33166 🌐 /mercoframes
✉ sales@mercoframes.net 🌐 www.mercoframes.com
☎ 305-882-0120 📱 Whatsapp 🛒 www.mercoframesusa.com