

Pre liew


LCD ACUITY SYSTEM














Randomization function helps prevent

Animation video for pediatric fixation.

Patient education images.

ETDRS, specialty testing charts.



As technology continues to influence your practice, it should equally influence your patients. The new PreView LCD Acuity System from Marco streamlines your refracting process while also impressing your patients. Equipped with all the computerized tests needed for a comprehensive eye exam, the new PreView delivers controlled results while embracing the most important feature—simplicity. Distinctively different and extraordinarily simple, experience the efficiency of computerized acuity testing with the new Preview from Marco.

PREVIEW LCD Acuity Device with Remote SPECIFICATIONS

ALL IN ONE COMPUTER

Type: 18.5 inch Monitor

Resolution: 1366 X 768 pixels (16:9 aspect ratio)

Brightness: 230 ±30 cd/m²

• CHARTS

VA Chart Type: Snellen Letters, Numbers, Children,

Tumbling E's

Masks: Single letter, Horizontal Line, Vertical Line

(Randomization)

Range: 20 feet (600 – 10)

Special Chart Type: R/G, Dots, ETDRS, Phoria with fixation,

Coincidences, Binocular, Stereo Worth 4 Dot Test, Astigmatic Clock Dial, Child Fixation Animation, Patient Education

Screens

Distance of Installation: 6 ft - 30 ft (limited upward optotype

after 21 ft.)

Chart Selection Device: Wireless IR remote control

Power supply:

Input requirements: 100-240V~ (1.5A) 50-60Hz Output:

19VDC 3.42A

Power consumption: 65VA

Dimensions and Weight: 18"(W) x 15"(H) x 1.5"(D)/10lbs.

Standard Accessories: Wall Mounting Plate and Fasteners















