


With Tuttnauer Kwiklave™, you only need one autoclave for all your sterilization needs.

From handpieces to large instrument cassettes, the Tuttnauer Kwiklave will sterilize in just 9 minutes.

- Preset cycles for unwrapped or wrapped instruments, liquids and drying time.
- Electronic control panel - displays temperature, sterilization time, drying time and pressure (PSI).
- One touch operation - after the start button is depressed, the electronics do the rest, including venting the steam away from any electronic components.
- Eliminates the need for specialized autoclave.
- Low water cut-off - prevents burnout of heating elements.
- Available in 7, 9, or 10 inch chambers, to suit all your needs.



Panel with optional printer.



Available as Models 2340EKP and Model 2540EKP.

MODEL	CHAMBER DIMENSIONS	OVERALL DIMENSIONS	TRAY DIMENSIONS	NO. OF TRAYS SUPPLIED	POWER CONSUMPTION	POWER REQUIREMENT	TOTAL CYCLE TIME	WEIGHT
MODEL 1730EK	7" dia x 12" D (170 dia x 300Dmm)	16" D x 17" W x 12.5" H (410D x 435W x 315Hmm)	4.8" W x 12" L x 1" H (120W x 300L x 25Hmm)	3	12 Amp. 1350W	110V 60Hz	From cold start 12 min. (Ster temp. 273°F/134°C) From Hot Start 9 Mins.	53 lbs./ 24kgs.
MODEL 2340EK Model 2340EKP equipped with optional printer	9" dia x 18" D (230 dia x 450Dmm)	20" D x 20" W x 15" H (500D x 495W x 375Hmm)	7" W x 16" L x 1" H (170W x 400L x 25Hmm)	3 (Holds 2 full size plus 2 half size IMS Cassettes)	10 Amp. 2200W	220V* 60Hz	From cold start 12 min. (Ster temp. 273°F/134°C) From Hot Start 9 Mins.	72 lbs./ 33 kgs.
MODEL 2540EK Model 2540EKP equipped with optional printer	10" dia x 18" D (250 dia x 450Dmm)	20" D x 20" W x 15" H (500D x 495W x 375Hmm)	7" W x 16" L x 1" H (170W x 400L x 25Hmm)	4 (Holds 3 full size plus 3 half size IMS Cassettes)	10 Amp. 2200W	220V* 60Hz	From cold start 12 min. (Ster temp. 273°F/134°C) From Hot Start 9 Mins.	76 lbs./ 35 kgs.

*Models 2340EK and 2540EK require a full 220 volts in order to achieve sterilization for unwrapped instruments in 9 minutes from a hot start. Should your voltage be lower, safe sterilization can be attained by increasing the sterilization time.

AUTHORIZED DEALER



MERCOFRAMES OPTICAL CORP

5555 Nw 74 Ave. Miami, Fl 33166 305-882-0120
www.mercoframes.com sales@mercoframes.net

tuttnauer™

WHEN IT COMES TO SPEED, CAPACITY AND VALUE, NOTHING BEATS A KWIKLAVE.™

UNWRAPPED

WRAPPED

LIQUIDS

CAPACITY

Speed

tuttnauer 2340EK

VALUE

If you find yourself racing through a busy day, you'll really appreciate the Tuttnauer™ Kwiklave. The sterilization chamber in the Kwiklave reaches the required temperature and pressure in just a fraction of the time. Due to Tuttnauer's state-of-the-art technology, the sterilization cycle is reduced to approximately nine minutes.

The Kwiklave offers the full capacity and the sterilization assurance of a conventional Tuttnauer autoclave. In comparison, other

short cycle autoclaves accommodate only a small number of instruments, while others are so technically complex, they constantly malfunction. But the genius of the Kwiklave is that it was designed simply to be the best. And, because it's a Tuttnauer, it already has a proven track record.

The Tuttnauer Kwiklave is available in various chamber sizes and with an optional printer. Call your Tuttnauer dealer for complete information.



Shown: Model 2340EKP
with optional printer

tuttnauer™



MERCOFRAMES OPTICAL CORP

5555 Nw 74 Ave. Miami, FL 33166 305-882-0120

www.mercoframes.com sales@mercoframes.net