



## Light Guide Tip Configurations




**Standard Light Guide, 20 ga**  
Reorder No. TR9800-11, TR9800-61




**45° Light Pick, 20 ga**  
Reorder No. TR9801-11, TR9801-61




**45° Retinal Knife, 20 ga**  
Reorder No. TR9806-11, TR9806-61



**Wide Angle Light Guide, 20 ga**  
Reorder No. TR9821-11, TR9821-61



**Siepser Shielded Probe, 20 ga**  
Reorder No. TR9879-11



**Iris Retractor, 20 ga**  
Reorder No. TR9809-11



**DCR Probe, 20 ga**  
Reorder No. TR9880-11